

Biomechanics of Stand-Alone Lateral Interbody Fixation at L4-5 and L5-S1*

Kingsley R. Chin, M.D., *Institute for Modern and Innovative Surgery, Fort Lauderdale, FL*
L. Perez-Orribo, *Barrow Neurological Institute, Phoenix, AZ*
Philip M. Reyes, BSE, *Barrow Neurological Institute, Phoenix, AZ*
Anna G.U. Sawa, M.S., *Barrow Neurological Institute, Phoenix, AZ*
Steven C. Anagnost, M.D., *The Orthopaedic Center, Tulsa, OK*
Vivek P. Kushwaha, M.D., *Houston Orthopedic and Spine Hospital, Houston, TX*
Josue P. Gabriel, M.D., *St. Anthony's Memorial Hospital, Effingham, IL*
S. Craig Meyer, M.D., *Columbia Orthopaedic Group, Columbia, MO*
Carl A.R. Bruce, M.D., *University Hospital of the West Indies, Kingston, Jamaica*
Warren D. Yu, M.D., *George Washington University Hospital, Washington DC*
Neil R. Crawford, Ph.D., *Barrow Neurological Institute, Phoenix, AZ*




Research Performed At:
Barrow Neurological Institute, Spinal Biomechanics, Phoenix, AZ 85013

Principal Investigator:
Neil R. Crawford, Ph.D.



Co-Investigator:
Kingsley R. Chin, M.D.
Institute for Modern and Innovative Surgery,
Fort Lauderdale, FL 33311



Background Context

Transection of lateral annulus and removal of intervertebral disc are destabilizing procedures. However, it is unclear to what extent the insertion of a lateral interbody cage into the disc space restores stability to the motion segment. This question is of interest at L4-L5 where this is currently being done clinically as well as at L5-S1, where techniques for lateral interbody placement are inevitably being developed.

Purpose


The goal of this *in vitro* study was to quantify the stability imparted by a stand-alone lateral interbody lumbar PEEK cage that spans the width of the vertebral body.

Study Design/Setting

Nondestructive *in vitro* flexibility test with paired comparison of stability in two conditions.

Patient Sample

Seven human cadaveric L3-S1 specimens were studied. Procedures were performed at L4-L5 and L5-S1.




S-LIFT (Sagittal Lateral Interbody Fusion Technology) (SpineFrontier, Inc.) MA, USA

Outcome Measures

The range of motion (ROM) was quantified at L4-L5 and L5-S1 during flexion, extension, axial rotation, and lateral bending.

Methods

Nonconstraining nondestructive pure moments (7.5 Nm maximum) were applied to specimens while recording angular L4-L5 and L5-S1 motion optoelectronically. Specimens were tested intact and then again after discectomies at L4-L5 and L5-S1 and insertion of lateral interbody spacers (S-LIFT) spanning the entire width of the endplate. Difference in mean ROM



in the two conditions at either level was assessed using paired two-tailed Student's t-tests.

Results

At L4-L5, insertion of S-LIFT resulted in mean reductions in ROM of 26% during flexion ($p=0.06$), 23% during extension ($p=0.09$), 15% during axial rotation ($p=0.19$) and 8% during lateral bending ($p=0.67$) relative to the intact condition. At L5-S1, insertion of S-LIFT resulted in mean reductions in ROM of 43% during flexion ($p=0.017$), 35% during extension ($p=0.016$), 16% during axial rotation ($p=0.30$) and 14% during lateral bending ($p=0.50$) relative to the intact condition.

Discussion & Conclusion

The surgeon can expect a stand-alone lateral interbody device to provide a reasonably stabilizing effect but not enough to impart rigidity. Therefore posterior supplemental fixation is needed for greater rigidity.

References

- Cappuccino A, Cornwall GB, Turner AW, Fogel GR, Duong HT, Kim KD, Brodke DS. Biomechanical Analysis and Review of Lateral Lumbar Fusion Constructs. *Spine*. 2010 Dec 15;35(26 Suppl):S361-7.
- Laws CJ, Coughlin DG, Lotz JC, Serhan HA, Hu SS. Direct Lateral Approach to Lumbar Fusion is a Biomechanically Equivalent Alternative to the Anterior Approach: An *In Vitro* Study. *Spine*. 2011 Sep 30.

*Pending Publication